



**Celtarys**  
RESEARCH



# **CELT-050 VHL E3 Ligase**

## **TR-FRET protocol**

# Protocol for E3 ubiquitin ligase VHL affinity binding assay using **TR-FRET**

CELT-050, a potent VHL E3 ligase fluorescent ligand, was used as fluorescent probe in a TR- FRET assay to study the affinity of a known inhibitor of E3 ubiquitin ligase VHL, VH298.

## Equipment, Materials and Reagents:

- Plate Reader with a TR-FRET optic module capable of emitting in the 390nm range and detecting in the 650nm range.
- TR-FRET compatible plates such as CulturPlate-384.
- Centrifuge with adaptor compatible with 384 plates.
- Recombinant His Tagged Human VBC complex.
- Anti-His Europium labeled Antibody.
- CELT-050, VHL E3 ligase fluorescent ligand

## Assay Buffer:

- 20 mM HEPES pH 7.5
- 150 mM NaCl
- 1 mM TCEP
- 0.01 % TWEEN20

## Master Mix:

In order to streamline the process and run as many compounds as possible, preparation of a master mix is recommended. This mixture should have 2 times the concentration of the final concentration in the well:

- 400 nM CELT-050.
- 100 nM His-VBC complex.
- 1 nM AntiHis-EuAb.
- Assay Buffer.

## Assay Procedure

The assay was performed using 384 Microplates with 15  $\mu$ L as a final assay volume.

The assay was performed using Europium (Anti-His-EuAb) as donor and His-VBC complex as protein of interest to be degraded.

A 2X compound concentration template plate is recommended diluting the compounds in the assay buffer. The final DMSO concentration in the wells should be 2 %.

7.5  $\mu$ L of each 2X template well are deposited in their corresponding wells, and once all compounds are ready the plate is centrifuged.

Afterwards, 7.5  $\mu$ L of the Master Mix are added to each well, totaling 15  $\mu$ L in each well. The plate is centrifuged again.

It is recommended to cover the plate from the light while incubating. Measures were taken every half an hour for 3 hours, as the signal stabilizes around the 90 minute mark.


## Assay Performance

·Z´ prime of 0.87.

·Different percentages of DMSO were tested, from 0.1 % to 2 % (final concentration in the wells). Up to 2 % is acceptable although the assay window is reduced with the increasing percentages of DMSO.

## Results

The figure below shows how increasing concentrations of VH298 compete with CELT-050 ligand and thereby prevent TR-FRET from occurring.



VH298 displayed the expected potency in good correlation with the literature (Frost et al., 2016, Potent and selective chemical probe of hypoxic signalling downstream of HIF- $\alpha$  hydroxylation via VHL inhibition. Nature Communications volume 7, Article number: 13312).