



## **CELT-503 General Protocol**

### **REAGENTS**

- **WASH BUFFER\*:**
  - Tris-HCl or HEPES (pH 7.4)
- **INCUBATION BUFFER\*:**
  - Tris-HCl or HEPES (pH 7.4), MgCl<sub>2</sub>, NaCl, GDP, Saponin
- **MEMBRANE (10x)\***
- **COMPOUNDS/CONTROL (10x)\***
- **LIGAND (10x): CELT-503 (Fluorescent GTP $\gamma$ S)**
  - Recommended Final Concentration: 100nM
    - Intermediate Aliquot: 1 $\mu$ M in Incubation Buffer

### **ASSAY CONDITIONS**

- Final Volume: 250 $\mu$ L/well
- Incubation\* in Incubation Plate with shaking at 25°C

### **ASSAY PROTOCOL**

1. Add 175 $\mu$ L/well of Incubation Buffer to the Incubation Plate
2. Add 25 $\mu$ L/well of Membrane
3. Add 25 $\mu$ L/well of Compounds/Control
4. Incubate for 15 minutes
5. Add 25 $\mu$ L/well of CELT-503
6. Incubate for 1 hour
7. Transfer 200 $\mu$ L to the Filter Plate
8. Filter
9. Wash 3 times with 200 $\mu$ L/well of Wash Buffer
10. Read at  $\lambda$ Ex: 589 –  $\lambda$ Em: 616

*\*Optimize concentrations and incubation conditions for each receptor (GPCR) for optimal results*

For more details on the setup of assays using GTP $\gamma$ S, you can consult the following protocol:

[https://resources.perkinelmer.com/corporate/content/relatedmaterials/posters/sps\\_trfbasedgtpbindingassay.pdf](https://resources.perkinelmer.com/corporate/content/relatedmaterials/posters/sps_trfbasedgtpbindingassay.pdf)